

Autoclave (with pump vacuum)

Before the sterilization cycle starts, a pump generates a powerful depression or vacuum, removing all the air pockets still present in the sterilization chamber, and inside the wrappings and materials themselves.

The vacuum pump is constructed with special membranes for the aspiration of hot air and steam. Pump-generated vacuum ensures optimum quality of the sterilization cycle for both wrapped and unwrapped instruments.

At the end of the sterilization cycle, the pump starts up again to facilitate and improve the drying phase.



Height: 385mm

Width: 450mm

Depth: 540mm

Weight: 50kg

2053-S

- 5 Operating cycles (134°C - 121°C);
- capacity to sterilize S-type loads;
- provision for connection to the water mains;
- provision for memory card connection and for PC connection;
- built-in printer.

